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6. Polishing bur, color-coded with a yellow ring

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6.2 Speed range

The polishing bur is used at a speed of 20,000 rpm. Please use slight pressure only; this will result in con-siderable increase of the service life of the bur and create a highly brilliant surface on all precious metal alloys.

6.3 Cleaning the polishing bur

Use a soft cloth or steam to clean the polishing bur.

Note: To ensure consistent quality of the results, damaged and bent instruments need to be replaced im-mediately. Damaged instruments increase the risk of injury.

7. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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Instructions for use

Tools for milling techniques

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Indication range

The rotary tools for milling techniques are used in all fields of dental-technical milling work.

2. Indication

Processing metal surfaces in the milling unit/system.

3. Hazard and safety information

When working with tools for milling techniques, safety goggles and standard safety clothing must be worn. The maximum speed for the tools, which is indicated on the packaging, must be observed and must not be exceeded.

4. Wax bur

4.1 Modelling the wax crown

All modelling waxes suitable crown and bridge restorations can be used for modelling.

4.2 Speed range

Use the wax bur within a range of 3,000 to 5,000 rpm. It is recommended to exert only slight pressure when milling wax.

4.3 Cleaning the wax bur

The especially developed cutting geometry of the wax bur produces razor-sharp edges. Therefore wax residues should not be removed with the fingers to avoid any injuries. Use a brush, soft cloth or steam to remove wax residues.

5. Profile bur, groove bur and shoulder bur

5.1 Milling and drilling oil

To achieve smoother running of the bur on the surface of the material, a smoother milled surface and pro-longed service life, it is urgently recommended to use milling and drilling oil, (REF 550 0000 8). Thanks to considerably reduced evaporation, the oil produces a homogeneous lubricating film on the milled surface. It is not recommended to use essential (etheral) oils.

5.2 Speed range

The profile bur, groove bur and shoulder bur are used for milling within a range of 15,000 to 20,000 rpm. Please use slight pressure only; this will result in considerable increase of the service life of the bur and create a slight brilliance on the surface.

5.3 Cleaning the profile bur, groove bur and shoulder bur

Use a soft cloth or steam to clean the profile bur, groove bur and shoulder bur.

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